



**MEDIT**  
**i500**  
intraoral scanner

The Medit i500 intraoral scanner is a win-win-win for you, your lab and your patient. It saves time, money, and impression materials. Your patients will be impressed with the i500's advanced technology as well as enjoy improved comfort over traditional impression methods.

With its impressive accuracy, speed and powderless system, the i500 allows for a smoother scanning experience, shortening turnaround times and increasing the productivity of your clinic.



**Spec Data**

**Small Tip & Lightweight**

Tip Size: 19 x 15.2 mm  
Hand piece Weight: 276g

**High Speed & Accuracy**

Accuracy: Single Die < 10 Micron  
Full Arch < 100 Micron  
Scan time: Quadrant = 2-3 Min  
Full Arch = 4-6 Min

**Comfortable & Productive**

Scanning Volume: 14 x 13 mm  
Scan Distance: Touch & Hover

**Open & Convenient**

Open System  
Powder Free

**HD Video & Real Color**

3D in motion video  
3D full color streaming

**Disinfection & Ortho readiness**

Tip 20x Autoclavable  
Ready for all Orthodontic uses

“It provides a scan experience that rivals some of the top scanners on the market and only for a fraction of the price.”

-Institute of Digital Dentistry


**Medit i500**

Cerec/Sirona Omnicam

Carestream CS3600

TRIOS 3 Wireless

**Scanner Comparison Chart**



Scanner Comparison Chart	Medit i500	Cerec/Sirona Omnicam	Carestream CS3600	TRIOS 3 Wireless
Tip Size	19x15.2mm	16x16.5mm	20x17mm	23x20mm
Wand Weight	276 g	311 g	326 g	340 g
Annual Fee	None	\$199/month	None	\$125-\$375/month
Scan Method	Triangulation	Triangulation	Triangulation	Confocal Laser
Full-Arch Scan Time	5 min	8-10 min	8-10 min	5 min
Cloud Services	Included	Optional*	N/A	Included*
Case Sharing	Open	Closed	Open	Simi-Open*
Workflow Management	Integrated	Optional*	No	Yes
MSRP	\$17,995	\$30,000	\$36,000	\$48,900

\*Additional fees may apply